

TELE PACK + COMPACT HYSTEROSCOPY




The All-in-One Solution for Uterine Diagnosis and Treatment

The TELE PACK + compact endoscopy solution is an investment that expands the depth and breadth of services you provide to your patients. Designed to meet the diverse diagnostic and treatment needs of your Gynecology office, this unique platform is compatible with the full range of cost-effective reusable hysteroscopic solutions from KARL STORZ.

- Compact, all-in-one solution fits even the smallest exam room
- Provides better patient care, efficiency, and increased satisfaction
- EHR integration makes it easy to communicate and track patient progress
- See it. Capture it. Treat it.

TELE PACK + COMPACT HYSTEROSCOPY

The All-in-One Solution for Uterine Diagnosis and Treatment



- + 18.5" touch-screen with on-screen keyboard
- + Full HD
- + Hygienic glass surface



- + LED light source
- + Automatic and manual light intensity modes
- + Stroboscopy



- + Multi-specialty for ENT, Urology, Gynecology, and more
- + Compatible with rigid, flexible, and single-use scopes



- + Integrated recording and playback
- + 50 GB internal storage
- + Encrypted data



- + Cyber secure
- + Network capability
- + User roles

Ordering Information	
TP101	TELE PACK + , endoscopic video unit, including Power Cord
TH110	IMAGE1 S™ HX One-Chip FULL HD Camera Head , high-definition
495NA	Fiber Optic Light Cable , 3.5 mm x 230 cm
TP002	Camera Head/Irrigation Bottle Holder , for use with TELE PACK +
20040282	USB Flash Drive , 32 GB
9700TPX	Rolling stand
Hysteroscope Ordering Information	
High-Flow Single Sheath – One hysteroscope that delivers superior flow rate, efficiency and comfort	
26153EA	High-Flow Stereoscope Sheath , size 5mm
26161VS	Examination Sheath , diameter 4.1 mm
26120BA	HOPKINS® Forward-Oblique Telescope 30° , diameter 2.9 mm, length 30 cm
26020FA	HOPKINS® Telescope 12° , diameter 2.9 mm, length 30 cm
CAMPO TROPHYscope® – Compact and versatile all-in-one solution	
26008BAC	CAMPO TROPHYscope® , HOPKINS® telescope 30°, size 2.9 mm, length 24 cm
26152DA	Continuous-Flow Examination Sheath , size 3.7 mm
26152DB	Continuous-Flow Operating Sheath , size 4.4 mm
BETTOCCHI® Systems (available in 4 mm, 5 mm, 6 mm, and Extended length versions) – ask for details	
5 Fr. Semi-rigid Mechanical Reusable Hand Instruments	
26159UHW/26160UHW	Biopsy and Grasping Forceps , semi-rigid, 5 Fr., 34 cm or 40 cm
26159EHW/26160EHW	Scissors , semi-rigid, blunt, 5 Fr., 34 cm or 40 cm
26159H/26160H	HESSELING Tenaculum Grasping Forceps , 5 Fr., 34 cm or 40 cm
26159BHW/26160BHW	Biopsy Spoon Forceps , semi-rigid, 5 Fr., 34 cm or 40 cm
Accessories	
26153EAU/10	Single-Use Sealing Cap – 10 per pack
260814FM-5	Fluid Management Pack – 5 pack/box

* Integrated recording and playback and networking in the U.S. will require a future software release, and all networking functionalities will also require the StreamConnect® media management application.

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com

STORZ
KARL STORZ – ENDOSKOPE
THE DIAMOND STANDARD

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com

Caution: Federal (USA) law restricts this device to sale by or on the order of physician.

© 2020 KARL STORZ Endoscopy-America, Inc.

97128009 GYN 1 1.0 us 06/2020/SE-E-US