

INTRAUTERINE BIGATTI SHAVER (IBS®)



Small Access, Great Effect

No thermal damage to tissue thanks to the purely mechanical removal of polyps, myomas or placental remnants









- Direct suction of resected tissue possible
- Effective, precise and safe alternative to resectoscopes
- Minimal outer diameter, optimized procedure time*, enhanced safety, accessibility and patient satisfaction*
- Fully reusable hysteroscopic shaver system
- NEW: More intraoperative possibilities thanks to various 4 mm instruments

* <http://www.fvvo.be/assets/725/FVVinObGyn-10-115.pdf>



INTRAUTERINE BIGATTI SHAVER (IBS®)

With an optical system that features an outer diameter of only 6 mm, the Intrauterine BIGATTI Shaver (IBS®) is inserted into the minimally dilated cervix. Thanks to the 4 mm shaver blade, the immediate removal of intrauterine pathologies is now possible outside the OR, e.g., in office settings and clinical outpatient departments or IVF centers.

Telescope and Instruments	
26208AMK	HOPKINS® Telescope 6° , Set
26702050	DRILLCUT-X® II Shaver Handpiece GYN
26208SA	Shaver Blade GYN , straight, sterilizable, concave cutting edge, double serrated, oval cutting window, Ø 4 mm
26208SB	Shaver Blade GYN , straight, sterilizable, double serrated cutting edge, rectangular cutting window, Ø 4 mm
26208SZ	Coagulation Electrode , bipolar, for use with Intrauterine BIGATTI Shaver (IBS®)
26176LE	Bipolar High Frequency Cord , length 300 cm

Further intraoperative options with various 4 mm instruments				
				
	26410MDM	26410MDN	26410MDO	CLICKLINE Dissecting and Grasping Forceps
	26410MVM	26410MVN	26410MVO	CLICKLINE Dissecting and Grasping Forceps
	26410MWM			CLICKLINE Scissors
	26410SAM			CLICKLINE Scissors STERILE 

Motor Unit	
26701001-1	UNIDRIVE® S III, Gyn Set , motor system, power supply 100-120/230-240 VAC, 50/60 Hz
20701070	Control cable

Continuous-Flow Pump System, comprising two ENDOMAT® SELECT pumps		
Suction	UP210	ENDOMAT® SELECT
	UP603	IBS® Shaver Software
	030647-10*	Tubing Set for Suction STERILE 
Irrigation	UP210	ENDOMAT® SELECT
	UP602	HYSTEROSCOPY Software
	UP610	Advanced Package Software
	031523-10*	Tubing Set for Irrigation, pressure-controlled STERILE 
Accessories	030970-10*	Tissue Trap Filter



For direct contact email us at gyn@karlstorz.com

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com

STORZ
KARL STORZ – ENDOSKOPE
THE DIAMOND STANDARD

www.karlstorz.com

