


# THE BIPOLAR 5 mm OFFICE RESECTOSCOPE



Click or scan here  
to watch the teaser video

## Versatile use









- Minimal outer diameter
  - Optimized intervention time\*
  - Increased safety, accessibility and patient satisfaction\*
  - Possible to perform resection in the doctor's office or clinical outpatient setting, even without anesthesia
  - Little or even no dilation required (depending on the patient)
- Optional use as „see & treat“ hysteroscope
- Real bipolar system with current not returned via the sheath




\* <http://www.fvvo.be/assets/725/FVWinObGyn-10-115.pdf>.


# THE BIPOLAR 5 mm OFFICE RESECTOSCOPE

Due to its small diameter, this resectoscope offers the possibility to perform resections outside the OR. Infertility patients or those suffering from abnormal bleeding can thus be treated in the doctor's office or clinical outpatient setting, even without anesthesia.

Optic and Working Element Set	
26120 AA	<b>HOPKINS® Straight Forward Telescope 0°</b> , diameter 2.9 mm, length 30 cm
26053 SCK	<b>Resectoscope Sheath</b> , 15 Fr., oblique beak, rotating inner sheath, color code: green
26053 OC	<b>Standard Obturator</b>
26053 EBH	<b>Working Element Set</b> , bipolar
26053 CD	<b>Telescope Bridge</b> , with channel for semi rigid 5 Fr. operating instruments, for use with Resectoscope Sheath 26053SC

HF Electrodes for various indications (sterile, for single use, package of 10)			
	011050-10*	<b>Cutting Loop</b> , bipolar, 15 Fr.,	STERILE 
	011051-10*	<b>Cutting Electrode</b> , bipolar, pointed, 15 Fr.,	STERILE 
	011052-10*	<b>Coagulation Electrode</b> , bipolar, ball end, 15 Fr.,	STERILE 
	011053-10*	<b>Cutting Loop</b> , bipolar, longitudinal, 15 Fr.,	STERILE 

Cold Loops for blunt dissection/enucleation or retrieving resected tissue		
	26053R	MAZZON <b>Cold Loop</b> , straight, 15 Fr.
	26053T	MAZZON <b>Cold Loop</b> , rake-shaped, 15 Fr.
	26053U	MAZZON <b>Cold Loop</b> , knife-shaped, 15 Fr.

HF Device		Pump System	
UH400	<b>AUTOCON® III 400</b> , high-end HF generator, without CE certification, power supply 220--240 VAC, 50/60 Hz	UP210	<b>ENDOMAT® SELECT</b>
UH400E	<b>Same</b> , with CE certification	UP602	<b>HYSTEROSCOPY Software</b>
		031523-10*	<b>Tubing Set for Irrigation</b> , pressure-controlled <span style="float: right;">STERILE </span>



For direct contact email us at [gyn@karlstorz.com](mailto:gyn@karlstorz.com)

It is recommended to check the suitability of the product for the intended procedure prior to use.  
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

KARL STORZ SE & Co. KG  
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
Postbox 230, 78503 Tuttlingen/Germany  
Phone: +49 7461 708-0  
Fax: +49 7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)

[www.karlstorz.com](http://www.karlstorz.com)

**STORZ**  
KARL STORZ — ENDOSKOPE  
THE DIAMOND STANDARD

