


KARL STORZ PEDIATRIC RIGID BRONCHOSCOPY CONTAINER SYSTEM



Protect Your Pediatric Bronchoscopy Instrument Investment

Our KSZ-39316F Container System is specially designed to protect your Pediatric Rigid Bronchoscopy Instrument Sets during transport, sterilization, and storage.



- Validated and approved for steam sterilization
- Save time and money over blue wrap
- Anodized to prevent pitting and corrosion





STORZ
KARL STORZ – ENDOSKOPE

KARL STORZ PEDIATRIC RIGID BRONCHOSCOPY CONTAINER SYSTEM

Protect Your Pediatric Bronchoscopy Instrument Investment

The KSZ-39316F Sterilization Container System is designed to house your KARL STORZ pediatric rigid bronchoscopy instruments through transport, sterilization and storage. It offers cost-effective, hassle-free sterilization for your valuable pediatric rigid bronchoscopy instrument investment. KARL STORZ sterilization systems are made from aircraft-grade anodized aluminum to prevent pitting and corrosion. They also help eliminate the need for blue sterilization wrap, while saving time and money.

Ordering Information	
 	<p>KSZ-39316F</p> <p>Container System, for pediatric bronchoscopy instruments. Includes two instrument racks, base basket, container base, and lid. 23.5" x 11" x 8".</p>

Accessories For Steam Sterilization	
	<p>KSZ-SCF01</p> <p>Singe-Use Paper Filters for Steam Sterilization, 1000 per case. 100% Cellulose and water repellent, 7.5" round.</p>
	<p>KSZ-SCL01</p> <p>Steam Sterilization Dual Process Indicator ID Cards, 1000 per case.</p>
	<p>KSZ-SCI001</p> <p>Steam Sterilization Dual Process Indicator Strips, 1000 per case. A dual purpose indicator strip with room for process identification. 8" long with perforation to reconfigure up to two 4" long indicator strips.</p>
	<p>KSZ-SCS01B</p> <p>Steam Sterilization Tamper Evident Seals, 1000 per pack.</p>

Contact your local KARL STORZ representative for more information about these or other KARL STORZ products, or call 800-421-0837.

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
 Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
 Postbox 230, 78503 Tuttlingen/Germany
 Phone: +49 7461 708-0
 Fax: +49 7461 708-105
 E-Mail: info@karlstorz.com
 www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
 2151 East Grand Avenue
 El Segundo, CA 90245-5017, USA
 Phone: +1 424 218-8100
 Phone toll free: 800 421-0837 (US only)
 Fax: +1 424 218-8525
 Fax toll free: 800 321-1304 (US only)
 E-Mail: communications@karlstorz.com