

# KARL STORZ DUAL YOKE HOSE — ELIMINATE THE RISK OF RUNNING ON EMPTY



## KARL STORZ Dual Yoke Hose

Designed to drive efficiency in the operating room, the KARL STORZ Dual Yoke Hose eliminates the distraction of exchanging CO<sup>2</sup> cylinders during your surgical procedure. This dual yoke assembly makes the process of changing to a full CO<sup>2</sup> cylinder 30 times faster\*. Avoid the disruption of running out of CO<sup>2</sup> and maintain continuous visualization of your surgical image with the KARL STORZ Dual Yoke Hose.

- Connects two CO<sup>2</sup> E-Cylinder Tanks to one insufflator
- Guarantees quick access to a backup CO<sup>2</sup> cylinder for extended laparoscopic procedures
- On and Off labels are carved into the system for quick setting identification
- Prevents loss of visualization due to loss of insufflation
- Backward compatible with the KARL STORZ 9606 series of surgical carts
- Never get caught without a backup

\*<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4432722/>

# KARL STORZ DUAL YOKE HOSE — ELIMINATE THE RISK OF RUNNING ON EMPTY


## KARL STORZ Dual Yoke Hose

Maintaining pneumoperitoneum throughout the surgical procedure is critical to sustain cavity visualization. With a traditional CO<sub>2</sub> cylinder system, you run the risk of running out of CO<sub>2</sub> and losing visualization. In fact, most surgeons have reported loss of surgical visualization during critical times in laparoscopic surgery due to the decrease in pneumoperitoneum during CO<sub>2</sub> exchange\*. To avoid these distractions, KARL STORZ has a dual yoke hose assembly for rapid swapping between CO<sub>2</sub> cylinders.

With just a twist of the cylinder selector switch, you can change from one empty cylinder while operating off a second CO<sub>2</sub> cylinder. Having the KARL STORZ Dual Yoke Hose as part of your visual and insufflation technology means you'll never get caught without a backup.

Ordering Information	
XUSM0007	Dual Yoke Assembly

Cart Compatibility Information	
COR Cart Series	UG130U-KT
9606 Series	9606-32
	9606A-32



\*<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4432722/>

**Contact your local KARL STORZ representative for more information  
about these or other KARL STORZ products, or call 800-421-0837.**

It is recommended to check the suitability of the product for the intended procedure prior to use.  
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

KARL STORZ SE & Co. KG  
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
Postbox 230, 78503 Tuttlingen/Germany  
Phone: +49 7461 708-0  
Fax: +49 7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
  
[www.karlstorz.com](http://www.karlstorz.com)



KARL STORZ Endoscopy-America, Inc.  
2151 East Grand Avenue  
El Segundo, CA 90245-5017, USA  
Phone: +1 424 218-8100  
Phone toll free: 800 421-0837 (US only)  
Fax: +1 424 218-8525  
Fax toll free: 800 321-1304 (US only)  
E-Mail: [communications@karlstorz.com](mailto:communications@karlstorz.com)