

KARL STORZ RotaTip
for bariatric surgery



RotaTip straight

KARL STORZ ROTATIP – for bariatric surgery



The ROTATIP line is a forceps and scissors system that can be dismantled into three parts. Thanks to the robust design of the outer sheath, utilization of this system can prove beneficial for bariatric surgery.



- Instruments feature an outer diameter of 5.5 mm
- Can be used with all standard 6 mm trocars
- Instruments have a length of 43 cm
- Instruments can be dismantled into three parts: Handle, outer sheath and working insert

4 jaw types are currently available:

for grasping



23410NR **ROTATIP Forceps Insert**, single action jaws, atraumatic, fenestrated



23410AFR **ROTATIP Forceps Insert**, double action jaws, atraumatic, fenestrated

for dissecting



23410MLR **ROTATIP Forceps Insert**, double action jaws, long

for cutting



23410MSR **ROTATIP Scissors Insert**, double action jaws, curved



The outer sheath features cm-markings (5 cm and 10 cm) and has been assigned the order number 23400.

- The outer diameter of this instrument is approx. 20% larger than outer sheaths in the CLICKLINE series. Consequently, this outer sheath is more robust and – in combination with the cm-marking – makes the ROTATIP system suitable for interventions on adipose or large patients.
- The handle has somewhat larger dimensions and the position of the rotary wheel has been optimized.
- No articulation of the joint mechanism when opening the jaws.

Outer Sheath, length 43 cm	Handle
23400	23151

Working Insert	Complete Instrument
Single action jaws	
23410ON	23451ONR
	ROTATIP Grasping Forceps , atraumatic, fenestrated
Double action jaws	
23410AF	23451AFR
	ROTATIP Grasping Forceps , atraumatic, fenestrated
23410ML	23451MLR
	ROTATIP Dissecting and Grasping Forceps , long
23410MS	23451MSR
	ROTATIP Scissors , curved

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.



75 Years

*Shaping the Future
of Endoscopy with you*



STORZ
KARL STORZ—ENDOSKOPE

THE DIAMOND STANDARD

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com

