

CLEANING BRUSHES FOR FLEXIBLE ENDOSCOPES



Don't Trust Your Endoscopes with Anything Less

Designed for effective and efficient cleaning, our single-use brushes ensure compliance with product IFUs and help maintain the integrity and precision of our flexible endoscopes.

- The only brushes designated for use with KARL STORZ products
- Gentle, thorough, confident cleaning
- Available in a wide range of shapes and sizes
- 10% discount on convenient bundled packages









CLEANING BRUSHES FOR FLEXIBLE ENDOSCOPES

Don't Trust Your Endoscopes with Anything Less

- KARL STORZ single-use cleaning brushes are specifically designed for use with KARL STORZ endoscopes—makes it easy to comply with our product IFUs
- The nonabrasive bristles provide a soft yet effective scrubbing action while protecting delicate instruments from damage
- The comprehensive range of styles and sizes covers our complete line of flexible endoscopes

Ordering is easy, with packages that include 50 of each of the brushes needed for a particular endoscope. Choose from six different bundles, and enjoy the convenience and savings that come with ordering our bundled packs at 10% off!

BRUSH-KT1	BRUSH-KT2	BRUSH-KT3	BRUSH-KT4	BRUSH-KT5	BRUSH-KT6
27652-05: 10 packs	27652-05: 10 packs	27652-05: 10 packs	27652-05: 10 packs	27652-05: 10 packs	27652-05: 10 packs
9990/50: 1 pack	9990/50: 1 pack	9990/50: 1 pack	9990/50: 1 pack	9990/50: 1 pack	9990/50: 1 pack
11276CL2/10: 5 packs	110910/50:1 pack	110910/50: 1 pack	110910/50: 1 pack	11276CL/10: 5 packs	110910/50: 1 pack
		11276CL2/10: 5 packs	11275CS/10: 5 packs	11275CS/10: 5 packs	110917-50: 1 pack

Featured Brushes	
	27652-05 Single use flat cleaning brush for the exterior cleaning. 5 per package. 15 mm width x 15.1 cm length x 7.5 mm height.
	9990/50 Single use cleaning brush for the suction wells. 50 per package. 11 mm in diameter x 5.5 cm in length.
	110910-50 Single use cleaning brush for working channels of 2.8 - 4.2 mm diameter. 50 per package. 5.0 mm in diameter x 220 cm length.
	110917-50 Single use cleaning brush for working channels of 1.2 mm. Double ended. 50 per package. 1.5 mm in diameter x 220 cm length.
	11275CS/10 Single use short cleaning brush. 10 per package. 4 mm in diameter and 10 cm in length.
	11276CL2/10 Single use long cleaning brush for working channels 1.2-1.5 mm. 10 per package. 2.5 mm in diameter and 110 cm in length.

Contact your KARL STORZ representative or call 1-800-421-0837 to place an order.

It is recommended to check the suitability of the product for the intended procedure prior to use. Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com