

KARL STORZ ADULT RIGID BRONCHOSCOPY CONTAINER SYSTEM



Protect Your Adult Rigid Bronchoscopy Instrument Investment

The KSZ-27717D Container System is specially designed to protect your Adult Rigid Bronchoscopy Instrument Sets during transport, reprocessing, and storage.


- Validated and approved for steam sterilization
- Save time and money over blue wrap
- Anodized to prevent pitting and corrosion





STORZ
KARL STORZ – ENDOSKOPE

KARL STORZ ADULT RIGID BRONCHOSCOPY CONTAINER SYSTEM

Protect Your Adult Rigid Bronchoscopy Instrument Investment

The KSZ-27717D Sterilization Container System is designed to house your KARL STORZ adult rigid bronchoscopy instruments through transport, sterilization and storage. It offers cost-effective, hassle-free sterilization for your valuable adult rigid bronchoscopy instrument investment. KARL STORZ sterilization systems are made from aircraft-grade anodized aluminum to prevent pitting and corrosion. They also help eliminate the need for blue sterilization wrap, while saving time and money.

Ordering Information		
	KSZ-27717D	Container System , for adult bronchoscopy instruments and telescopes. Includes base basket, accessory basket, two racks, container base, and lid. 27" x 11" x 8"

Accessories For Steam Sterilization		Description
	KSZ-SCF01	Disposable Paper Filters for Steam Sterilization , 1000 per case. 100% Cellulose and water repellent, 7.5" round.
	KSZ-SCL01	Steam Sterilization Dual Process Indicator ID Cards , 1000 per case.
	KSZ-SCI001	Steam Sterilization Dual Process Indicator Strips , 1000 per case. A dual purpose indicator strip with room for process identification. 8" long with perforation to reconfigure up to two 4" long indicator strips.
	KSZ-SCS01B	Steam Sterilization Tamper Evident Seals , 1000 per pack

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com



FOR U.S. USE ONLY

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com

© 2021 KARL STORZ Endoscopy-America, Inc.

Caution: Federal (USA) law restricts this device to sale by or on the order of physician.