


C-VIEW® DIGITAL CYSTOSCOPE



Cystoscopy Simplified

The new C-VIEW® has been simplified to provide the perfect balance of image quality and ease-of-use to meet the needs of high-volume Urology Offices.

Simple to Use:

- Plug-and-Go Usability
- Navigate Challenging Anatomy with Ease
- Designed for Durability

Simple to Maintain:

- Efficient Reprocessing
- Next-Day Repair/Exchange

C-VIEW® DIGITAL CYSTOSCOPE

Cystoscopy Simplified

Product Specifications

Outside Diameter	15.6 Fr.
Distal Tip	10.9 Fr.
Working Channel	7 Fr.
Deflection	210°/140°
Direction/Angle of View	0°/100°
Light Source	Integrated LED
Working Length	37 cm

*Contoured tip for smooth access
and increased patient comfort*



C-VIEW® Digital Cystoscope	
11272VUEK	C-VIEW® Digital Cystoscope , US deflection
11272VUEK-R	C-VIEW® Digital Cystoscope , reverse deflection
Flexible Instruments & Accessories	
27425FG	Grasping Forceps , F/Stone (stent removal), 5 Fr., 60 cm
27425ZG	Biopsy Forceps , flexible, double, 5 Fr., 60 cm
27723T	Bugbee Electrode , monopolar, ball, 4 Fr., 73 cm
27502	Luer Connector w/ Stopcock
27550ZAG	Self-Sealing Cap (single-use)
C-MOS Office/Mobile Video System	
8403ZX	C-MAC® Monitor/CCU
20290301	C-HUB® II CCU
VTG0011-UG	Sony 21" HD Monitor , with table-top stand
TELE PACK+ (All-In-One Endoscopic Imaging Platform)	
TP101	TELE PACK+ endoscopic video unit , including power cord
TH110 or TH111	IMAGE1 S™ HX Camera Heads , high-definition
9700TPX	Rolling Stand

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com