


CLICKLINE PRECISION INSTRUMENTS FOR COLORECTAL SURGERY



Hand Instruments Tailored to Your Needs

For a colorectal surgeon, the right instruments are crucial. Whether working to manipulate the bowel or remove a section of tissue, appropriate grasping and dissecting tools are essential. KARL STORZ is proud to offer an excellent line of MIS instruments and associated accessories:








- Bowel graspers in plain, Hunter, and Prestige styles
- Atraumatic graspers and dissectors in Universal, Bullet Nose, Wavy, and Maryland styles
- Teethed graspers with Tiger, Claw, and Alligator jaws
- Anvil graspers

CLICKLINE PRECISION INSTRUMENTS FOR COLORECTAL SURGERY





Hand Instruments Tailored to Your Needs

Here are the most commonly used forceps for colorectal procedures in standard size. Please see the KARL STORZ CLICKline brochure for more options, including different lengths and sizes.











Bowel Grasping Options

	33310C Insert for bowel grasper , fenestrated, double-action, 5 mm x 36 cm
	33310CC Insert for CROCE-OLMI grasping forceps , atraumatic, fenestrated, curved, single-action, 5 mm x 36 cm
	33310DY Insert for DeBAKEY grasping forceps , atraumatic, single-action, 5 mm x 36 cm
	33310DYD Same , double-action <i>*Similar to a Hunter Grasper</i>
	33310DK Insert for DeBAKEY grasping forceps , curved and slender jaws, atraumatic, double-action, 5 mm x 36 cm
	33310LF Insert for grasping forceps , atraumatic, hollow jaws, single action, 5 mm x 36 cm
	33310LFD Same , double action <i>*Similar to a Prestige Grasper</i>


Teethed Grasping Options

	33310MG Insert for grasping forceps , 2 x 4 teeth, for grasping and removal of solid organs, in particular during adhesiolysis, single-action, 5 mm x 36 cm
	33310G Insert for grasping forceps , 2 x 4 teeth, double-action, 5 mm x 36 cm <i>(Also known as Kocher or Tiger jaws in OR)</i>
	33310FM Insert for grasping forceps , 2 x 3 teeth, single-action, 5 mm x 36 cm <i>(Also known as Claw in OR)</i>
	33310HM Insert for dissecting and grasping forceps , double-action, 5 mm x 36 cm <i>(Also known as Alligator or Crocodile jaws)</i>

Atraumatic Grasping and Dissecting Options

	33310ON Insert for grasping forceps, with especially fine atraumatic serration , fenestrated, single-action, 5 mm x 36 cm
	33310OD Same , double-action <i>(Also known as Johan or Universal in OR)</i>
	33310UL Insert for REDDICK-OLSEN dissecting and grasping forceps , robust, double-action, 5 mm x 36 cm <i>(Also known as "Bullet" Nose in OR)</i>
	33310AF Insert for grasping forceps , atraumatic, fenestrated, double-action, 5 mm x 36 cm
	33310SN Insert for SCHNEIDER Lymph Node grasping forceps , atraumatic, single-action, 5 mm x 36 cm
	33310WT Insert for grasping forceps , atraumatic, single-action, 5 mm x 36 cm
	33310WTD Same , double-action <i>(Also known as Wavy grasper in OR)</i>
	33310MD Insert for KELLY dissecting and grasping forceps , double-action, 5 mm x 36 cm <i>(Also known as Maryland in OR)</i>
	33310ML Insert for dissecting and grasping forceps , long, double-action, 5 mm x 36 cm <i>(Also known as Maryland in OR)</i>
	33310R Insert for dissecting and grasping forceps , right-angled, double-action, 5 mm x 36 cm <i>(Also known as Mixer in OR)</i>

Anvil Graspers

	33310AV Insert for grasping forceps , for stapler pressure plates, double-action, 5 mm x 36 cm
--	---

It is recommended to check the suitability of the product for the intended procedure prior to use. Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com