


CLICKLINE PRECISION INSTRUMENTS FOR GYN LAP SURGERY



Hand Instruments Tailored to Your Needs







Gynecological laparoscopy requires a broad range of instruments. Whether working to preserve tissue or manipulate the hardest fibroids, surgeons need the right tools for the job. KARL STORZ is proud to offer an excellent line of MIS instruments and associated accessories:




- Atraumatic graspers and dissectors i.e. Maryland, Universal, and Bullet Nose styles
- Bowel graspers
- Tooth graspers with different modalities
- GYN-specific graspers i.e. Tenaculum, Biopsy Punch, and Oviduct Grasper





CLICKLINE PRECISION INSTRUMENTS FOR GYN LAP SURGERY






Hand Instruments Tailored to Your Needs

Here are the most commonly used forceps for GYN LAP procedures in standard size. Please see the KARL STORZ CLICKline brochure for more options, including different lengths and sizes.

Atraumatic Grasping and Dissecting Options	
	33310MD Insert for KELLY dissecting and grasping forceps , double-action, 5 mm x 36 cm <i>(Also known as Maryland in OR)</i>
	33310ML Insert for KELLY dissecting and grasping forceps , long, double-action, 5 mm x 36 cm <i>(Also known as Maryland in OR)</i>
	33310ON Insert for grasping forceps, with especially fine atraumatic serration , fenestrated, single-action, 5 mm x 36 cm
	33310OD Same , double-action <i>(Also known as Johan or Universal in OR)</i>
	33310UL Insert for REDDICK-OLSEN dissecting and grasping forceps , robust, double-action, 5 mm x 36 cm <i>(Also known as Bullet "Nose" in OR)</i>
	33310AF Insert for grasping forceps , atraumatic, fenestrated, double-action, 5 mm x 36 cm

Bowel Grasping Options	
	33310C Insert for bowel grasper , fenestrated, double-action, 5 mm x 36 cm
	33310LF Insert for grasping forceps , atraumatic, hollow jaws, single-action, 5 mm x 36 cm
	33310LFD Same , double-action <i>(Similar to a Prestige Grasper)</i>

Teethed Grasping Options	
	33310MG Insert for grasping forceps , 2 x 4 teeth for grasping and removal of solid organs, in particular during adhesiolysis, single-action, 5 mm x 36 cm
	33310G Insert for grasping forceps , 2 x 4 teeth, double-action, 5 mm x 36 cm <i>(Also known as Kocher or Tiger jaws in OR)</i>
	33310ME Insert for grasping forceps, with multiple teeth , width of jaws 4.8 mm, for atraumatic and accurate grasping
	33310FM Insert for claw forceps , 2 x 3 teeth, single-action, 5 mm x 36 cm

GYN-Specific Grasping Options	
	33310UM Insert for ovary grasping forceps , serrated, double-action, 5 mm x 36 cm <i>(Similar to a Tenaculum Grasper)</i>
	34310DZ Insert for biopsy punch forceps , with two teeth, single-action, 5 mm x 36 cm
	34310DZL Insert for biopsy punch forceps , with two teeth, long, single-action, 5 mm x 36 cm
	33310VA Insert for mesoovarian forceps , forceps, double-action, 5 mm x 36 cm
	33310VT Insert for oviduct forceps , double-action, 5 mm x 36 cm

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com

STORZ
KARL STORZ — ENDOSKOPE
THE DIAMOND STANDARD

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com