

CLICKLINE PRECISION INSTRUMENTS FOR GENERAL SURGERY



Hand Instruments Tailored to Your Needs












KARL STORZ is proud to offer an excellent line of MIS instruments and associated accessories. These graspers/dissectors are:





- Easy to disassemble and reassemble
- Less cost of repair due to interchangeable parts
- Modular that enables thorough cleaning and sterilization




CLICKLINE PRECISION INSTRUMENTS FOR GENERAL SURGERY

Hand Instruments Tailored to Your Needs

Here are the most commonly used forceps for LAP procedures in standard size procedures. Please see the KARL STORZ CLICKline brochure for more options, including different lengths and sizes.

Atraumatic Grasping and Dissecting Options	
	33310FG Insert for grasping forceps , atraumatic, double-action, 5 mm x 36 cm
	33310MD Insert for KELLY dissecting and grasping forceps , double-action, 5 mm x 36 cm <i>(Also known as Maryland in OR)</i>
	33310ML Insert for KELLY dissecting and grasping forceps , long, double-action, 5 mm x 36 cm <i>(Also known as Maryland in OR)</i>
	33310R Insert for dissecting and grasping forceps , right-angled, double-action, 5 mm x 36 cm <i>(Also known as Mixer in OR)</i>
	33310ON Insert for grasping forceps, with especially fine atraumatic serration , fenestrated, single-action, 5 mm x 36 cm
	33310OD Same, double-action jaws <i>(Also known as Johan or Universal in OR)</i>
	33310UL Insert for REDDICK-OLSEN dissecting and grasping forceps , robust, double-action, 5 mm x 36 cm <i>(Also known as "Bullet" Nose in OR)</i>
	33310LF Insert for grasping forceps , atraumatic, hollow jaws, single-action 5 mm x 36 cm
	33310LFD Same , double-action jaws <i>*Similar to a Prestige Grasper</i>
	33310DY Insert for DeBAKEY grasping forceps , atraumatic, single-action, 5 mm x 36 cm
	33310DYD Same , double-action jaws <i>*Similar to a Hunter Grasper</i>

Teethed Grasping Options	
	33310MG Insert for grasping and removal of solid organs , in particular during adhesiolysis, single-action, 5 mm x 36 cm <i>(Also known as Kocher or Tiger jaws in OR)</i>
	33310G Insert for grasping forceps , 2 x 4 teeth, double-action, 5 mm x 36 cm <i>(Also known as Kocher or Tiger jaws in OR)</i>
	33310MA Insert for grasping forceps , 1 x 2 teeth, single-action, 5 mm x 36 cm <i>(Also known as Cobra jaws in OR)</i>
	33310HM Insert for dissecting and grasping forceps , double-action, 5 mm x 36 cm <i>(Also known as Alligator or Crocodile jaws)</i>

Scissors Options	
	34310EH Insert for hook scissors , 5 mm x 36 cm
	34310MS Insert for METZENBAUM scissors , curved, 5 mm x 36 cm
	34310MT Insert for scissors , serrated, straight, 5 mm x 36 cm

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com

STORZ
KARL STORZ — ENDOSKOPE
THE DIAMOND STANDARD

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com