

HD VIDEO MEDIASTINOSCOPY FROM KARL STORZ



Systematic Simplicity

Video Mediastinoscopy - Plug & Play

In thoracic surgery, video mediastinoscopy is of great value for the histological clarification of mediastinal lymph nodes and tumor staging. Tissue structures are clearly displayed for the purpose of diagnosis and aid with precise dissection of lymph nodes.

Equipped with future-oriented CMOS technology and an LED light source integrated in the handle, the HD video mediastinoscope provides the best possible imaging using the modular KARL STORZ IMAGE1 S™ camera system and matching instruments. THINK ONE!





STORZ
KARL STORZ – ENDOSKOPE

HD VIDEO MEDIASTINOSCOPY FROM KARL STORZ

Future oriented – Compact – Straightforward

- Fully integrated, Single component mediastinoscope
- Improved Visualization for lymph node biopsy
- Reduced Set up time
- Compact and Autoclavable



Ordering Information		
	10973 HD	Video Mediastinoscope , with continuous proximal lateral slit, length 15 cm, for use with IMAGE1 S™
	10971 EN	Coagulation Suction Cannula , insulated, with connector pin for unipolar coagulation, angled handle, diameter 5 mm, length 20 cm
	33151 MK	CLICKLINE METZENBAUM Scissors , rotating, dismantling, insulated, with connector pin for unipolar coagulation, with Luer-Lock irrigation connector for cleaning, double action jaws, curved, size 5 mm, length 25 cm
	10971 GO	MEDIAFIT Biopsy Forceps , oval jaws, size 8 x 16 mm, sheath diameter 5 mm, length 21 cm

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com



FOR U.S. USE ONLY

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com

© 2021 KARL STORZ Endoscopy-America, Inc.

Caution: Federal (USA) law restricts this device to sale by or on the order of physician.