

WALL PLATE, 2N1 IP ENCODER

MULTI-FORMAT 4K IP ENCODER (WUIS3403)

SPECIFICATIONS
SHEET



HIGHLIGHTS

- A multi-format wall plate that converts 4K DisplayPort or HDMI 2.0 video signals into raw KARL STORZ NEO IP network packets. The NEO IP system is capable of distributing ultra high definition 4K video in pristine quality with imperceptible latency to the end user, which is critical when performing endoscopic surgeries.
- Designed in a compact/durable 1-gang wall plate. Ideal for wall and equipment boom locations in operating rooms
- Supports resolutions up to 4096x2160 @ 60Hz video sources
- Auto-sensing input video detection

FEATURE	DETAIL
Input Connections	DisplayPort 1.2, HDMI 2.0
Output Connections	Dual LC Fiber, 10Gb Network
Power	12VDC, 2A Use only KARL STORZ approved medical grade power sources.
Thermal	82 BTU / Hr.
Operating Temperature	32 - 104°F (0 - 40°C)
Operating Humidity	0-95% Non-Condensing
Dimensions Ext. (WxHxD)	2.9" x 4.7" x 3.6"
Applied Standards	IEC 60601-1, IEC 60601-1-2, AAMI ES60601-1, IEC 60825-1, RoHS3, REACH, FDA Class I Medical Devices

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

STORZ
KARL STORZ—ENDOSKOPE


FOR U.S. USE ONLY

WALL PLATE, 2N1 IP ENCODER

MULTI-FORMAT 4K IP ENCODER (WUIS3403)

SPECIFICATIONS
SHEET

INSTALLATION



NOTES

- Stub Conduit 2" min. and 4" MAX above trim ring
- The drawings above are generic and for coordination only, not for construction. The KARL STORZ Customer Success Manager assigned to your project should be consulted for any site specific requirements.

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com