

# Unidrive<sup>®</sup>



# Select

Power Anywhere



**STORZ**  
KARL STORZ — ENDOSKOPE

Discover the highlights of the new UNIDRIVE® SELECT motor system for ENT, neurosurgery and spine surgery:



**Power, precision and versatility combined in one system:**

✓ **Multidisciplinary and modular concept**

Various software modules for interdisciplinary use and customizable expansion of application fields

✓ **Coordinated system**

The new KS HIVE™ communication standard enables irrigation to be automatically initiated when the UNIDRIVE® SELECT is active. This provides a system that features an inter-coordinated pump and motor unit.

✓ **Design**

New KARL STORZ unit design with glass front and touch screen

✓ **One-handed control**

User-friendly tubing set thanks to easy insertion

✓ **Control**

New wireless LED footswitch for flexible positioning

Unidrive®  
 Select  
Power Anywhere



✓ **Broad motor and handpiece portfolio**

Handpieces for diverse fields of application in ENT, neurosurgery and spine surgery

✓ **Power**

Powerful motors/handpieces with up to 80,000 rpm

✓ **Integrated irrigation**

Some handpieces with integrated irrigation

✓ **Length-adjustable burrs**


The new high-speed burrs are adjustable to three lengths

✓ **New micro ear burrs**

Slim and light micro ear burrs for access to the middle ear

✓ **Shaver handpieces**

Compatibility with existing DRILLCUT-X® II shaver handpieces



It is recommended to check the suitability of the product for the intended procedure prior to use.  
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

More than  
**75**  
Years

*Shaping the Future  
of Endoscopy with you*

**STORZ**  
**KARL STORZ—ENDOSKOPE**

*THE DIAMOND STANDARD*

KARL STORZ SE & Co. KG  
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
Postbox 230, 78503 Tuttlingen/Germany  
Phone: +49 7461 708-0  
Fax: +49 7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
[www.karlstorz.com](http://www.karlstorz.com)

