

Permanent Availability

SCENARA® Content Management



Data availability with SCENARA® Content Management independent of time or location

SCENARA® Content Management is used to manage, edit and save the content generated along the entire patient pathway. This maintains an overview of images, videos and reports assigned to each patient. This information can be retrieved from any PC, tablet or smartphone at all times. Data security is guaranteed via encryption according to the latest security standards. A backup strategy prevents data loss and the user administration ensures that data can only be viewed by authorized personnel.

Image and video management

SCENARA® .MEDIA allows straightforward, intuitive processing and storage of endoscopic image and video data generated in the OR. The export function in various file formats such as PDF or PowerPoint provides support for the processing of a wide range of applications and is an optimal addition to OR documentation.

- Central access to image, video and patient data
- Direct image and video upload from a smartphone
- Comparison of one or more interventions via parallel viewing of two videos
- Adaptive streaming to ensure continuous streaming
- Rapid navigation in videos via a surgical timeline
- Ease of use thanks to fast forward and frame-by-frame navigation
- Video editing, snapshots and image processing with a few clicks
- Intuitive addition of annotations
- Data transfer to PACS





Procedure management

SCENARA® .PROCEDURES supports the planning of procedures and a smooth workflow in the OR. Through an interface with the hospital information system, patient and order data can be synchronized at any time. Changes to the workflow can be easily and quickly made and are immediately visible in all applications.

- Planning of the OR workflow
- Filter options and search menu
- HL7 interface
- DICOM Worklist Provider



Documentation

SCENARA® .REPORTS supports OR documentation through an efficient generation of standardized OR reports. Predefined protocols simplify orientation and ensure completeness of the data. Along the entire patient pathway, documented information and OR images can be entered manually or directly imported from **SCENARA® .TASKS** and **SCENARA® .MEDIA**, thus guaranteeing the best possible flexibility of use.

- Supports documentation of medical findings
- Linking of image and diagnostic data
- Easy PDF generation
- Data export to the hospital information system

More than
75
Years

*Shaping the Future
of Endoscopy with you*

STORZ
KARL STORZ—ENDOSKOPE

THE DIAMOND STANDARD

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com

