

C-MAC FIVE S 5.3 SUPERIOR IMAGE & HANDLING IN A SINGLE USE SCOPE

C-MAC[®]
CHOICE CHANGES EVERYTHING



97078044 AN 1 2.0 us 12/2021/SE-E-US

Choice Changes Everything

The new FIVE S 5.3 scope combines the best of both worlds: the convenience of single-use with the performance of reusables. This ground-breaking device delivers the consistency, simplicity, ergonomics, and image quality that you have come to expect from KARL STORZ.

Manufactured with demanding precision, the single-use FIVE S 5.3 achieves minimal variation where it matters most.

The FIVE S 5.3 fits seamlessly into the existing C-MAC[®] family and is compatible with the C-MAC[®], C-HUB, as well as TELE PACK+. The FIVE S 5.3 single-use scope adds one more device to the most comprehensive platform with the most comprehensive connectivity.

NEW order number: 0915612-06 (6 endoscopes per package)			
Distal tip outer diameter	5.3 mm	Deflection up/down:	180°/180°
Working channel diameter:	2.2 mm	Angle of view:	90°
Working length:	65 cm	Direction of view:	0°

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

KARL STORZ Endoscopy-America, Inc. 2151 East Grand Avenue El Segundo, CA 90245-5017/USA
www.karlstorz.com

FOR U.S. USE ONLY

STORZ
KARL STORZ – ENDOSKOPE