

AIDA® XPAND INTERFACE

SPECIFICATIONS SHEET

(WUIS3523, WUIS3524)

AIDA® Xpand Interface
(WUIS3523, WUIS3524)

AIDA® Image Capture
(WD300)

AIDA® SmartScreen®
(WM200)



FEATURES

- A single clutter-free interface that connects the AIDA WD300 image-capture and documentation system to the NEO IP operating room integration solution.
- 2nd generation video-over-IP chipset reduces power consumption, heat load and noise from cooling fans.
- Models for interfacing with single and dual channel image-capture systems.

FEATURE	WUIS3523 Single Channel	WUIS3524 Dual Channel
Power	70W	
Weight	9.6 LBS	
Dimensions (WxHxD)	Standalone: 12 × 2.7 × 13.9 in. / 305 × 66 × 353 mm - 6 in. Minimum Cable Clearance Behind Unit Mounted with WD300 AIDA and WM200 SmartScreen: 12 × 7.5 × 14.9 in. / 305 × 191 × 379 mm (as shown above)	
Thermal	239 BTU/HR	
Operating Temperature	32–104°F (0–40°C)	
Operating Humidity	0–95% Non-Condensing	
Compatibility	NEO IP with software v2.2.1.3 or later	
Applied Standards	IEC 60601-1, RoHS3, REACH, FDA Class I Medical Devices	
Included Accessories	3 SFP+ Modules, 1 Medical-Grade Power Cable	5 SFP+ Modules, 1 Medical-Grade Power Cable
Install Kit	WUIS3625 - Single Channel Install Kit for AIDA WD300-KT or WD350-KT Image-Capture System	WUIS3626 - Dual Channel Install Kit for AIDA WD300-KT or WD350-KT Image-Capture System

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue | El Segundo, CA 90245-5017, USA
Phone: +1 424-218-8100 (Toll free - US only:) 800-421-0837
Fax: +1 424-218-8525 (Toll free - US only:) 800-321-1304
E-Mail: communications@karlstorz.com

www.karlstorz.com

FOLLOW KARL STORZ

