

# COMPACT ENDOSCOPY CART & STAND SOLUTIONS DESIGNED FOR OFFICES AND CLINICS



## Endoscopy Cart & TELE PACK Stand Solutions for a Range of Procedural Applications

### **Maximize Space and Efficiency**

Both the [Compact Endoscopy Cart](#) and the [Compact Endoscopy Stand](#) feature an “open concept” minimal footprint designed to improve workflow in smaller treatment rooms at offices and clinics.

### **Adaptable for Various Treatment Options**



These holding solutions can be paired with a range of optional accessories, including an IV pole, reusable and single-use scope holders, and MORE to enable many configuration options so you can better serve your patients.

# COMPACT ENDOSCOPY CART & STAND SOLUTIONS DESIGNED FOR OFFICES AND CLINICS

## Endoscopy Cart & TELE PACK Stand Solutions for a Range of Procedural Applications

**NEW** Compact Endoscopy Mounting Solutions designed to meet the endoscopic video storage needs for an array of procedural applications in the office and clinic setting.

- Compact footprint to maximize space and efficiency in smaller treatment rooms
- Accessories enable adaptability for various treatment options
- Open design improves ventilation and communication
- Adjustable display height for an improved viewing experience
- Covered cable management channel allows for a clear appearance

Compact Endoscopy Cart & Stand	Part No.
 <p><b>Compact Endoscopy Cart</b> includes:</p> <ul style="list-style-type: none"> <li>• Keyboard Tray (sliding)</li> <li>• 5.5" Drawer</li> <li>• Camera Holder</li> <li>• 8 Outlet Isolation Transformer</li> <li>• Cable Management Wraps</li> </ul>	XUSM0001
 <p><b>Compact Endoscopy Stand</b> includes:</p> <ul style="list-style-type: none"> <li>• Keyboard Tray</li> <li>• Basket</li> <li>• Cable Management Wraps</li> </ul>	XUSM0013

Optional Accessories	Part No.
IV Pole	XUSM0001-4
Videoscope Holder	XUSM0001-5
Single Use Scope Holder Basket	XUSM0001-21
Proshield Tube Holder	XUSM0001-23
Camera Holder	XUSM0001-6
Suction Canister Holder	XUSM0001-7
Dual CCD Scope Holder*	XUSM0001-18
Adjustable Shelf	XUSM0001-1

\* Dual CCD Scope Holder is intended to hold CCD Flexible GI Silver Scopes®, and is not designed to hold CMOS videoscopes. This accessory also requires the IV Pole accessory, P/N XUSM0001-4 to be mounted.

**Contact your local KARL STORZ representative for more information about this or other KARL STORZ products and programs, or call 800-421-0837.**

**FOR U.S. USE ONLY**

It is recommended to check the suitability of the product for the intended procedure prior to use.  
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.



**KARL STORZ Endoscopy-America, Inc.**  
2151 East Grand Avenue | El Segundo, CA 90245-5017, USA  
Phone: +1 424-218-8100 (Toll free - US only:) 800-421-0837  
Fax: +1 424-218-8525 (Toll free - US only:) 800-321-1304  
E-Mail: [communications@karlstorz.com](mailto:communications@karlstorz.com)

[www.karlstorz.com](http://www.karlstorz.com)

FOLLOW KARL STORZ

