

ENDOMAT® SELECT

UP220

The Choice is Yours



Discover the **Highlights** of the **ENDOMAT® SELECT:**



✓ **Flexibility thanks to modularity**

Various software modules for interdisciplinary use and customizable expansion of application fields


✓ **Design**

New KARL STORZ unit design with glass front and touch screen

✓ **One-handed control**

User-friendly application thanks to easy-to-insert cassette tubing set

✓ **Intuitive user interface with height difference functionality**



Independent, flexible and free positioning on the tower possible thanks to height difference functionality

✓ **Coordinated with the UNIDRIVE® SELECT system**

Thanks to the new KS HIVE® communication standard, the irrigation cycle of the ENDOMAT® SELECT is automatically initiated when the UNIDRIVE® SELECT is active. This provides a system that features an inter-coordinated system pump and motor unit for ENT, neurosurgery and spine surgery.



New Software License UP611

New Tubing Set 031531-10



It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

More than
75
Years

*Shaping the Future
of Endoscopy with you*