

INTRAUTERINE BIGATTI SHAVER (IBS®)



Small Access, Great Effect

No thermal damage to tissue thanks to the purely mechanical removal of polyps, myomas or placental remnants. Resected tissue is directly aspirated and collected in a tissue trap filter.

- Effective, precise and safe alternative to resectoscopes*
- Minimal outer diameter for shorter polypectomy times compared to conventional resection via resectoscope*
- Fully reusable hysteroscopic shaver system
- More intraoperative possibilities thanks to various 4 mm instruments

* <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6658200/>

INTRAUTERINE BIGATTI SHAVER (IBS®)

The Intrauterine BIGATTI Shaver (IBS®) features an optical system with an outer diameter of just 6 mm and a conical obturator tip to facilitate insertion into the minimally dilated cervix.

Telescope and Instruments	
26208AMK	HOPKINS® Telescope 6° , Set
26702050	DRILLCUT-X® II Shaver Handpiece GYN
26208SA	Shaver Blade GYN , straight, sterilizable, concave cutting edge, double serrated, oval cutting window, Ø 4 mm
26208SB	Shaver Blade GYN , straight, sterilizable, double serrated cutting edge, rectangular cutting window, Ø 4 mm
26208SZ	Coagulation Electrode , bipolar, for use with Intrauterine BIGATTI Shaver (IBS®)
26176LE	Bipolar High Frequency Cord , length 300 cm

Further intraoperative options with various 4 mm instruments				
	26410MDM	26410MDN	26410MDO	CLICKLINE Dissecting and Grasping Forceps
	26410MVM	26410MVN	26410MVO	CLICKLINE Dissecting and Grasping Forceps
	26410MWM			CLICKLINE Scissors
	26410SAM			CLICKLINE Scissors

Motor Unit	
26701001-1	UNIDRIVE® S III, Gyn Set , motor system, power supply 100-120/230-240 VAC, 50/60 Hz
20701070	Control cable

Pump System, comprising two ENDOMAT® SELECT pumps		
Suction	UP210	ENDOMAT® SELECT
	UP603	IBS® Shaver Software
	030647-10	Tubing Set for Suction STERILE
Irrigation	UP210	ENDOMAT® SELECT
	UP602	HYSTEROSCOPY Software
	UP610	Advanced Package Software
	031523-10	Tubing Set for Irrigation, pressure-controlled STERILE
Accessories	030970-10	Tissue Trap Filter

** The ENDOMAT® SELECT is only intended to be used for suction in context of the IBS® Shaver indication. Hysteroscopy with continuous flow, i.e., pressure- and flow-controlled irrigation with simultaneous suction, requires the HAMOU® ENDOMAT® (item no. 26331101-1).

For direct contact email us at gyn@karlstorz.com

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.