

## SINGLE-USE VIDEO URETERO-RENSCOPE

### HIGHLIGHTS

- Single-use
- 105° Angle of view
- 0° Direction of view
- 270° Deflection (up and down)
- 70 cm Working length
- 3.5 Fr Working channel (3 o'clock position)
- 9 Fr Sheath size
- CMOS Chip
- Integrated LED light source



Ready When You Are



### SPECIFICATIONS

#### PRODUCT DETAIL

091271-06	Video Uretero-Renoscope FLEX-X <sup>C1</sup> , with <b>positive</b> deflection mechanism, package of 6 pcs
091279-06	Video Uretero-Renoscope FLEX-X <sup>C1</sup> , with <b>contra-positive</b> deflection mechanism, package of 6 pcs

#### IMAGE1 S™ COMPATIBILITY

TC028	E-BOX Videoendoscope Adaptor for X-LINK
TC200EN or TC201EN	IMAGE1 S CONNECT®*
TC301EN	IMAGE1 S X-LINK
TM343	32" 4K Monitor

\* FLEX-X<sup>C1</sup> requires the IMAGE1 S™ software version 4.3 or higher. For the latest software version please contact Technical Customer Service ([technicalsupport@karlstorz.com](mailto:technicalsupport@karlstorz.com)).

Other imaging compatibilities like TELE PACK+, TELECAM C3 and C-MAC® HD Monitor are available for FLEX-X<sup>C1</sup>. Please contact your KARL STORZ representative for more information.



Please scan the QR code to sign-in for a demo and to see FLEX-X<sup>C1</sup> in action

It is recommended that the suitability of the products for the intended procedure be checked prior to use. Please note that the products listed here may not yet be available in all countries due to differences in approval requirements.