

# Configuration Planning

## ENDOMAT® SELECT

The Choice is Yours





## SURGERY



For flow-controlled applications such as, for example, laparoscopy, thoracoscopy, proctology:

031524-10 **Tubing Set**, irrigation, FC, STERILE   
sterile, for single use, package of 10,  
for use with KARL STORZ HAMOU® ENDOMAT®  
26331120-1, ENDOMAT® SELECT UP210  
and ENDOMAT® SELECT UP220



Procedures that use a flow-controlled tubing set bear the adjacent mark.



For hysteroscopic tissue ablation with the IBS® shaver system or the removal of tissue chips during resection (TUR-B/P):

030647-10 **Tubing Set**, suction, STERILE   
DS, sterile, for single use, package of 10,  
for use with KARL STORZ ENDOMAT® SELECT  
UP210 and ENDOMAT® SELECT UP220



Procedures that use a tubing set for direct suction bear the adjacent mark.



For use in combination with the CALCUSON lithotripsy unit in urology:

031647-10 **Tubing Set**, suction, STERILE   
via bottle, sterile, for single use, package of 10,  
for use with KARL STORZ ENDOMAT® SELECT  
UP210 and ENDOMAT® SELECT UP220



Procedures requiring a suction bottle bear the adjacent mark.



For use with the CLEARVISION® functionality in ENT:

031529-10 **Tubing Set**, irrigation, STERILE   
CV, sterile, for single use,  
package of 10, for use with  
KARL STORZ ENDOMAT® SELECT UP210




Procedures that use a tubing set for the CLEARVISION® functionality bear the adjacent mark.

In the daily surgical routine in the OR, the ENDOMAT® SELECT delivers a powerful yet low-noise irrigation of fluids during interventions related to the following specialties:

- Laparoscopy
- Thoracoscopy
- Proctology

### Software:

- UP601 **SURGERY Software**, license, allows selection of the procedures “LAP”, “THOR” and “PROCTO”
- Promotion Pack (031524-10)** STERILE   
**Tubing Set**, irrigation, FC, sterile,  
for single use, package of 10



## (OFFICE) HYSTEROSCOPY

The ENDOMAT® SELECT provides pressure-controlled dilation of the uterine cavity for diagnostic and operative hysteroscopy in the office setting or the clinical outpatient department. Alternatively, the user can select suction-controlled irrigation.

### Software:


UP602 **HYSTEROSCOPY Software**, license, allows selection of the procedure "HYS"

Please select the required Promotion Pack: STERILE 

**Promotion Pack (031523-10) Tubing Set**, irrigation, PC, sterile, for single use, package of 10

or

**Promotion Pack (030647-10) Tubing Set**, suction, DS, sterile, for single use, package of 10




## Hysteroscopic Tissue Ablation with the Intrauterine BIGATTI Shaver IBS®

In conjunction with the IBS® shaver, the ENDOMAT® SELECT can be used for the controlled suction of resected polyps, myomas or placental remnants.

### Software:

UP603 **IBS® SHAVER**, software, license, allows selection of the procedure "IBS®"

**Promotion Pack (030647-10)**, STERILE   
**Tubing Set**, suction, DS, sterile, for single use, package of 10



### Required accessories:

20701070 **Control Cable**, connectors 2 x LEMO 5-pin 0° connectors, length 100 cm, for transmission of footswitch control between the UNIDRIVE® S III 20701020-1 and ENDOMAT® SELECT UP210 or ENDOMAT® SELECT UP220

## UROLOGY

The ENDOMAT® SELECT accommodates the wide range of applications in the field of urology, allowing seamless integration into existing systems and providing pressure-regulated irrigation for diagnostic and operative interventions.

### Software:

UP604 **UROLOGY Software**, license, allows selection of the procedures "CYST", "RES", "URS", "CALCUSON", and "PCN"

Please select the required Promotion Pack: STERILE 


**Promotion Pack (031523-10) Tubing Set**, irrigation, PC, sterile, for single use, package of 10

or

**Promotion Pack (030647-10) Tubing Set**, suction, DS, sterile, for single use, package of 10

or

**Promotion Pack (031647-10) Tubing Set**, suction, BS, sterile, for single use, package of 10



### Further accessories:

20014130 **One-Pedal Footswitch**, digital, one-stage


20701070 **Control Cable**, connectors 2x LEMO 5-pin 0° connectors, length 100 cm, for transmission of footswitch control signal between the ENDOMAT® SELECT UP210 and UNIDRIVE® S III 20701020 such as between the ENDOMAT® SELECT UP210 and CALCUSON 27610020 or ENDOMAT® SELECT UP220 and CALCUSON 27610020

## ARTHROSCOPY

In the field of arthroscopy, the ENDOMAT® SELECT regulates pressure-controlled irrigation of fluids to maintain joint distension. Alternatively, the user can select suction-controlled irrigation.

### Software:

- UP605 **ARTHROSCOPY Software**, license, allows selection of the procedures "KNEE", "HIP", "SHOULDER" and "SMALL JOINTS", for use with the ENDOMAT® SELECT UP210
- Please select the required Promotion Pack:** STERILE 
- Promotion Pack (031523-10) Tubing Set**, irrigation, PC, sterile, for single use, package of 10
- or
- Promotion Pack (030647-10) Tubing Set**, suction, DS, sterile, for single use, package of 10



### Further accessories:

- 20014130 **One-Pedal Footswitch**, digital, one-stage
- UP006 **Control Cable**, connectors 1x LEMO 5-pin 0° and 1x LEMO 12-pin 0° connector, length 100 cm, for transmission of control signal between the UNIDRIVE® S III ARTHRO SCB 28723020-1 and the ENDOMAT® SELECT UP210 or the ENDOMAT® SELECT UP220

## ENT


The ENDOMAT® SELECT includes the following software module for irrigation of the lens.

### Software:

- UP606 **ENT Software**, allows selection of the procedure "CLEARVISION®"
- Promotion Pack (031529-10)** STERILE 
- Tubing Set**, irrigation, CV, sterile, for single use, package of 10, for use with the KARL STORZ ENDOMAT® SELECT UP210

### Required accessories:

- 20014230 **One-Pedal Footswitch**, digital, two-stage




## SPINE SURGERY

The ENDOMAT® SELECT offers the following software module for the irrigation of fluids during interventions in spine surgery.

### Software:

- |                                     |                                   |   |
|-------------------------------------|-----------------------------------|---|
| <input type="checkbox"/>            | UP607                             | <b>SPINE Software</b> , license, allows selection of the procedures “LUMBAR” and “THORACIC” |
| <input checked="" type="checkbox"/> | <b>Promotion Pack (031523-10)</b> | <b>Tubing Set</b> , irrigation, PC, sterile, for single use, package of 10                  |



## ADVANCED

In order to activate extended settings and additional functions, the “Advanced Package” can be selected as an additional\* software package.

### Software:

- |                          |       |   |
|--------------------------|-------|---|
| <input type="checkbox"/> | UP610 | <b>Advanced Package</b> , software, license |
|--------------------------|-------|---|

\* The advanced package only functions when at least one other software package is available.

The following features are available with the installation of the ADVANCED Software UP610:

- Adjustment of pressure (100/300/500 mmHg) during SURG applications
- Adjustment of inflow up to 3500 ml/min during SURG applications
- Adjustment of inflow (200/400/600 ml/min) during HYS, SPINE & URO applications
- Adjustment of inflow (1500/2000/2500 ml/min) during ARTHRO applications
- Possible to adjust the duration of the BOOST mode during ARTHRO & URO applications
- Adjustment of height to improve pump performance
- Choice of pressure unit, mmHg or cm H<sub>2</sub>O
- SCB integration

It is recommended to check the suitability of the product for the intended procedure prior to use.  
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

More than  
**75**  
Years

*Shaping the Future  
of Endoscopy with you*