

FLEX-X^{C1} SINGLE-USE URETEROSCOPE

SPECIFICATIONS SHEET

HIGHLIGHTS

- CMOS Chip
- Integrated LED light source
- 105° Angle of view
- 0° Direction of view
- 270° Deflection (up and down)
- 700 mm Working length
- 3.5 Fr. Inner channel
- 3 O'clock position working channel
- 9 Fr. Sheath size
- 9 Fr. Tapered distal tip



PRODUCT DETAIL

091279-06	Video Ureteroscope, sterile single use. Contrapositive (US) deflection, package of 6
091271-06	Video Ureteroscope, sterile single use. Positive (Reverse) deflection, package of 6

VIDEO COMPATIBILITY*

IMAGE1 S™ TC028	E-BOX (IMAGE1 S™ Single-Use Videoendoscope Adaptor)
TC200US or TC201US	CONNECT
TC301US	X-LINK
TM343	Monitor

*Other options available. Contact your KARL STORZ representative for more information.



It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue | El Segundo, CA 90245-5017, USA
Phone: 424-218-8100 Toll-free: 800-421-0837
Fax: 424-218-8525 Toll-free: 800-321-1304
E-Mail: communications@karlstorz.com

www.karlstorz.com

FOLLOW KARL STORZ



Caution: Federal (USA) law restricts this device to sale by or on the order of physician.

97138028 URO 1 1.0 US 01/2022/TS-E-US