

VITOM® 3D – MICROMANIPULATOR INTERFACE FOR CO₂ LASER MICROSURGERY




The VITOM® 3D delivers an ergonomic visualization experience supporting precise surgery, communication, and team effectiveness.


The MICROMANIPULATOR INTERFACE connects the VITOM® 3D to the CO₂ Laser Micromanipulator, enabling an ergonomic, upright sitting position and ensuring easy maneuverability of equipment.

STORZ
KARL STORZ – ENDOSKOPE

VITOM® 3D FOR CO₂ LASER MICROSURGERY



Compact design creates more work space for surgeon and staff.



The VITOM® 3D is the only exoscope on the market that enables an application with free beam lasers. The MICROMANIPULATOR INTERFACE is compatible with the Micromanipulators for microspot-based transoral surgery.











Integrates seamlessly with the IMAGE1 S™ modular scalable video system.



Improves OR staff engagement with a surgical display shared with the entire room.

Ordering Information

	TH004	Micromanipulator Interface VITOM® 3D, with clamping cylinder for mounting CO ₂ laser micromanipulators, for use with the VITOM® 3D
	TH200	VITOM® 3D, with zoom and focus function, integrated illumination and camera orientation correction, working distance 20–50 cm, fiber optic light transmission incorporated, wipe-disinfectable, for use with with the IMAGE1 S D3-LINK® (TC302) and the IMAGE1 PILOT (TC014)
	TC014	IMAGE1 PILOT, control unit with 3D wheel, 4 programmable function keys and USB port, for intuitive control of camera systems and connected units, for use with IMAGE1 S™ and the VITOM® 3D TH200
	UR500B/C UR002	ARTip CRUISE®, motorized holding arm, for short-term flexible positioning and dynamic movement of compatible KARL STORZ units during open, microsurgical and minimally invasive procedures VITOM® 3D Adaptor for the ARTip CRUISE®
	8668A	HINNI Distending Operating Laryngoscope, with right and left side wings, length 15 cm, for use with Smoke Evacuation Tube 8668M
	8890A	DEDO Operating Laryngoscope, with extra light handle, length 18 cm
	8790A	HAVAS Operating Laryngoscope, with extra light handle, large model, length 17 cm, with connector for Light Carrier 497HC
	8588A	WEERDA Distending Operating Laryngoscope, length 18 cm, large adult size

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.