

AUTOCON® III 300



The Standard Solution for Multidisciplinary High Frequency Surgery

The AUTOCON® III 300 HF generator is the latest member of the KARL STORZ product family for HF surgery. With a maximum output of unipolar 400 Watt and bipolar 200 Watt, it offers all the basic functions required in electrosurgery. One unipolar and one bipolar socket are available for electrosurgical cutting and coagulation in various specialties.

- High-resolution touch display
- Interdisciplinary use thanks to presettings and individually adjustable modes for various specialties
- Suitable for applications in arthroscopy
- Can be combined with the S-PILOT® for automatic smoke evacuation
- Can be combined with the AUTOCON® III ARGON for argon-assisted surgery

STORZ
KARL STORZ – ENDOSKOPE

AUTOCON® III 300

UH300	AUTOCON® III 300 power supply 220-240 VAC, 50/60 Hz, including mains cord
UH300E	Same, with CE certification
UH300U	AUTOCON® III 300 power supply 100-127 VAC, 50/60 Hz, including mains cord
UH300UE	Same, with CE certification



Accessories

UF901	One-Pedal Footswitch, with button for switchover function, for use with HF generators
UF902	Two-Pedal Footswitch, with button for switchover function, for use with HF generators
26005M	Unipolar High Frequency Cord, 300 cm (listed as example)
26176LE	Bipolar High Frequency Cord, 300 cm (listed as example)
27802	Neutral Electrode, for single use, contact surface divided into two, A=169 cm ² , package of 50
27806US	Neutral Electrode Connecting Cable, for use with Neutral Electrode 27802, length 500 cm
26520043	Electrode Handle, with two buttons for activating the unipolar generator (High Frequency Cable 26520045 required)
26520045	High Frequency Cable, for Electrode Handle 26520043, length 500 cm
UG910	COR HF Mobile Stand, optional accessories: UG810 Top Cover, UG815 Equipment Rail, UG902 Two-Pedal HF Footswitch Holder

Compatible equipment

UH500	AUTOCON® III ARGON, for use with AUTOCON® III 400 and AUTOCON® III 300 for non-contact coagulation and cutting
UH501	S-PILOT®, the solution for smoke evacuation, regulated suction and filtering of surgical smoke during endoscopic procedures, automatic activation with AUTOCON® III 400 and AUTOCON® III 300

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.