



HIGHLIGHTS

SCENARA® .doorsign offers detailed and automatically updated live information about the room and procedure status. It allows an overview of activities in the operating room without disrupting the OR workflow and features additional LED visualization of the OR status on the doorsign. Furthermore, the easy-to-use touchscreen provides a view of the entire daily schedule of the OR.

In addition to SCENARA® .doorsign it is possible to use SCENARA® .dashboard for a more enhanced overview of the current status of all ORs on a special touchscreen in a hospital corridor or at any other required location.

SPECIFICATIONS

PANEL	
LCD Size	10.1 inch
Resolution	1280 x 800 IPS
Brightness	500 cd/m2
Contrast Ratio	600:1
Viewing Angle	Horizontal 160°/Vertical 160°
Touch Type	Projected Capacitive (5 points)
Display Aspect Ratio	16:10

HIGHLIGHTS

SPECIFICATIONS

SYSTEM	
CPU	Cortex A17, Quad, Core, RK3288
GPU	MALI T764, 3D Graphic Engine
System Memory	2GB DDR3 SDRAM
Internal Storage	16GB eMMC Flash
Networking	LAN RJ45/WiFi/BT version 4.2
PoE Support	IEEE 802.3af and 802.3at (advised) PoE switch

PHYSICAL	
Housing	Plastic Housing
Color	White
Weight	610 g
Dimension	255 x 181 x 25 mm
Mounting	VESA 75
Power Supply	12V/2A
Power Usage	10 Watt, Maximum 19.5 Watt

FEATURES	
Surround LED Bar for Room Signage	
Power over Ethernet (PoE) and WLAN	
True Flat Projected Capacitive Multi-Touch	
HD (1280 x 800) IPS Panel	
Cortex A17, Quad Core, up to 2.0 GHz	
Commercial Grade (24/7 usage)	

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.