

KARL STORZ SENATOR Vacuum Pump

A successful, cost-efficient solution in a new design

The new KARL STORZ SENATOR vacuum pump is suitable for suctioning fluid or gas during endoscopic applications in doctors' offices, clinics, and hospitals.



It is easy to change locations thanks to the mobile suction solution or the integrated carrying handle so that the system can be used across specialties in various rooms.

In addition to easy handling, which plays a major role in the everyday OR routine, an uncomplicated combination with existing products is important to us. In addition to a stand-alone solution, the SENATOR vacuum pump can also be used in conjunction with the s-pilot® – our solution for smoke evacuation management – or the ENDOMAT® SELECT pump. This combination completes the imaging chain and ensures optimal viewing conditions in endoscopy.

STORZ
KARL STORZ – ENDOSKOPE

KARL STORZ SENATOR Vacuum Pump

A successful, cost-efficient solution in a new design

- High suction capacity of 30 l/min
- Hydrophobic bacterial filter
- Intuitive use
- Low noise and low vibration levels
- Can be used in various disciplines
- Easy-to-read manometer
- Rapid vacuum build-up
- Easy-to-clean surface
- For mounting to standard rails

Order numbers:



UP031 KARL STORZ SENATOR Vacuum Pump, Set

consisting of:

UP030

SENATOR Vacuum Pump, incl. mains cord,
1 connecting tube, 3 bacterial filters

UP032

Secretion Application Set, 2 l

consisting of:

UP030-5000168

Secretion Bottle, 2 l

UP030-500008

Bottle Cap for Secretion Bottle

UP030-500002

Patient Tube, 1.7 m, unsterile,
connection between secretion bottle and
instrument, for use with UP030, UP030C,
25320020 and 25320020C



Accessories

UP030-500001

Connecting Tube, short, 0.5 m, unsterile,
connection between filter and secretion bottle, for use with
UP030, UP030C, 25320020 and 25320020C

UP030-500084

Float for Overflow Protection

UP030-500002

Patient Tube, 1.7 m, unsterile,
connection between secretion bottle and instrument for use with
UP030, UP030C, 25320020 and 25320020C

UP030-5000168

Secretion Bottle, 2 l

UP030-500005

Hydrophobic Bacterial Filters

UP030-500008

Bottle Cap for Secretion Bottle

KARL STORZ SENATOR Vacuum Pump

A successful, cost-efficient solution in a new design

The KARL STORZ SENATOR vacuum pump can also be offered as a mobile suction unit on a small Mobile Stand UG911 from our proven COR series. The corresponding containers can be mounted on the stand and thus transported between OR rooms or for disposal.

Order numbers:



UG911	COR Pump Mobile Stand consisting of:
UG803	Subrack COR Mobile Stand, low
UG816	Pump Holder for COR Mobile Stand, for use with vacuum pumps KARL STORZ SENATOR and UNIMAT® 30

Accessories for KARL STORZ SENATOR

Setup with 1 Suction Liner

030350-05	Clamp Holder, standard
030301-01	Secretion Container, 1.5 l
030302-01	Secretion Container, 2.5 l
030315-30	Suction Liner, 1.5 l, with solidifier
030325-30	Suction Liner, 2.5 l, with solidifier
030310-40	Suction Liner, 1.5 l
030320-40	Suction Liner, 2.5 l
030362-10	Coupling Piece
030249-10	Tubing Set for Suction, patient side, 250 cm, with silicone tubing, sterile, for single use, package of 10
030248-10	Tubing Set for Suction, patient side, 250 cm, with male LUER-Lock, sterile, for single use, package of 10
030748-01	Tubing Set for Suction, unit side, 150 cm, unsterile

Option 1 Additional Components for Setup with 2 Suction Liners

030331-100 / 030332-01	Serial Connection, Patient, connection to suction liner (before filter), for single use
030341-01	Serial Connection, Vacuum, connection to suction liner (behind filter)
030351-01	Changeover Valve, clamp holder and tubing set

Option 2 Additional Components for Setup with Several (up to 5) Suction Liners

030331-100 / 030332-01	Serial Connection, Patient, connection to suction liner (before filter), for single use
030341-01	Serial Connection, Vacuum, connection to suction liner (behind filter)
030380-01	Trolley

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.



KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34 | 78532 Tuttlingen/Germany
Postbox 230 | 78503 Tuttlingen/Germany
Phone: +49 7461 708-0 | Fax: +49 7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com

FOLLOW KARL STORZ

