













ENT INSTRUMENT SET




Hemostasis



Photo	Product Number	Description/Specifications	Length	Size
	841111	Bipolar Forceps , insulated, straight, blunt, for use with Bipolar High-Frequency Cords 847000V	11 cm	
	841116	Bipolar Forceps , insulated, straight, pointed, for use with Bipolar High-Frequency Cords 847000V	16 cm	
	841216	Bipolar Forceps , insulated, straight, blunt, for use with Bipolar High-Frequency Cords 847000V	16 cm	1 mm
	841219	Bipolar Forceps , insulated, straight, blunt, for use with Bipolar High-Frequency Cords 847000V	19 cm	1 mm
	842016	Bipolar Forceps , insulated, tips curved, very delicate tips, for use with Bipolar High-Frequency Cords 847000V	16 cm	0.5 mm
	842116	Bipolar Forceps , insulated, tips curved, for use with Bipolar High-Frequency Cords 847000V	16 cm	0.5 mm
	842216	Bipolar Forceps , insulated, tips curved, blunt, for use with Bipolar High-Frequency Cords 847000V	16 cm	1 mm
	842319	Bipolar Forceps , insulated, tips curved, blunt, for use with Bipolar High-Frequency Cords 847000V	19 cm	1 mm
	842219	Bipolar Forceps , insulated, tips curved, blunt, for use with Bipolar High-Frequency Cords 847000V	19 cm	2 mm
	843016F	Bipolar Forceps , for part not insulated, tips curved, for use with Bipolar High-Frequency Cords 847000V	16 cm	0.2 mm
	843016	Bipolar Forceps , insulated, angular, very delicate, for use with Bipolar High-Frequency Cords 847000V	16 cm	0.2 mm
	843116	Bipolar Forceps , insulated, angled, pointed, for use with Bipolar High-Frequency Cords 847000V	16 cm	
	843219	Bipolar Forceps , insulated, angular, blunt, for nose, with integrated suction channel for epistaxis, with cut-off hole, for use with Bipolar High-Frequency Cords 847000V	19 cm	
	844116	Bipolar Forceps , insulated, bayonet shape, pointed, or use with Bipolar High-Frequency Cords 847000V	16 cm	
	844216	Bipolar Forceps , insulated, bayonet shape, blunt, for use with Bipolar High-Frequency Cords 847000V	16 cm	1 mm
	844219	Bipolar Forceps , insulated, bayonet shape, blunt, for use with Bipolar High-Frequency Cords 847000V	19 cm	1 mm
	844319	Bipolar Forceps , insulated, bayonet shape, slender, for use with Bipolar High-Frequency Cords 847000V	19 cm	0.5 mm
	844320	Bipolar Forceps , insulated, with integrated suction channel and cut-off hole, bayonet, blunt, for use with Bipolar High-Frequency Cords 847000V	20 cm	1.4 mm
	844321	Bipolar Forceps , insulated, with integrated suction channel and cut-off hole, bayonet, blunt, for use with Bipolar High-Frequency Cords 847000V	20 cm	2 mm

Photo	Product Number	Description/Specifications	Length	Size
	844423	Bipolar Forceps , insulated, bayonet, blunt, for use with Bipolar High-Frequency Cords 847000V	23 cm	0.3 mm
	844523	Bipolar Forceps , insulated, bayonet, blunt, for use with Bipolar High-Frequency Cords 847000V	23 cm	0.7 mm
	844623	Bipolar Forceps , insulated, bayonet, blunt, for use with Bipolar High-Frequency Cords 847000V	23 cm	1.2 mm
	845219	Bipolar Forceps , insulated, bayonet shape, tips curved, blunt, for use with Bipolar High-Frequency Cords 847000V	19 cm	2 mm
	845518	Bipolar Coagulation Forceps , insulated, for nasopharynx, tip curved 40°, blunt, for use with Bipolar High-Frequency Cords 847000V	18 cm	2 mm
	845721	Bipolar Coagulation Forceps , insulated, tip long curved 90° for nasopharynx and adenoidectomy, with integrated suction channel, with cut-off hole, for use with Bipolar High-Frequency Cords 847000V	21 cm	
	845821	Bipolar Coagulation Forceps , insulated, tip short curved 45° for oral cavity and tonsillectomy, with integrated suction channel, with cut-off hole, for use with Bipolar High-Frequency Cords 847000V	21 cm	
	839330	BRINER Bipolar Suction Cannula , angular, insulated, length of electrodes 3.5 mm, with cut-off hole, for use with Bipolar High-Frequency Cords 847000V	11 cm	4.5 mm
	847000V	Bipolar High-Frequency Cord , for Valleylab & Bovie Generators	300 cm	
	461010	STAMMBERGER Bipolar Suction Forceps , 15° upturned, for bipolar coagulation in paranasal areas, with suction channel, for use with Bipolar High-Frequency Cord 847002V	12.5 cm	
	461015	STAMMBERGER Bipolar Suction Forceps , 45° upturned, for bipolar coagulation in paranasal areas, with suction channel, for use with Bipolar High-Frequency Cord 847002V	12.5 cm	
	847002V	Bipolar High-Frequency Cord , for Valleylab coagulator, with two 2 mm cable sockets for Bipolar Suction Forceps 461010, 461015 and Bipolar Forceps 8615 A/AS	450 cm	

Clip Applicators (for use in Nose)

Photo	Product Number	Description/Specifications	Length
	452650A	MONTGOMERY-YOUNGS RHINOFORCE® Clip Applicator , for endonasal endoscopic sphenopalatine artery ligation, with suction channel, handle with spring, straight, with cleaning connector, for use with Ethicon Titanium Clips LT200	13 cm
	452650C	MONTGOMERY-YOUNGS RHINOFORCE® Clip Applicator , for endonasal endoscopic sphenopalatine artery ligation, with suction channel, handle with spring, jaws angled to the right, with cleaning connector, for use with Ethicon Titanium Clips LT200	13 cm
	452650D	MONTGOMERY-YOUNGS RHINOFORCE® Clip Applicator , for endonasal endoscopic sphenopalatine artery ligation, with suction channel, handle with spring, jaws angled to the left, with cleaning connector, for use with Ethicon Titanium Clips LT200	13 cm

TAKE-APART® Bipolar Forceps (for use in Nose & Neuro)



TAKE-APART® Bipolar Forceps, configuration consisting of:
 1x 26284HM-NEURO: Handle
 1x 26284AS-HN: Outer Sheath
 1x 26284BS-HN: Inner Sheath
 1 pack sterile, single-use inserts, available in various configurations below:




Photo	Product Number	Description/Specifications	Length	Size
	26176LV-HN	Bipolar High-Frequency Cord for Valleylab , reusable	300 cm	
	26176VLD-10-HN	Bipolar High-Frequency Cord for Valleylab , sterile, for single use, package of 10	300 cm	
	28164MGD-10	Bipolar High-Frequency Cord for MALIS , sterile, for single use, package of 10	300 cm	
	26284AS-HN	TAKE-APART® Bipolar Forceps outer sheath for 28164 series forceps	20 cm	3.4 mm
	26284BS-HN	TAKE-APART® Bipolar Forceps outer sheath for 28164 series forceps	20 cm	3.4 mm
	26284HM-NEURO	Bipolar Ring Handle		

Hemostasis (for use in Larynx)






Photo	Product Number	Description/Specifications	Length	Size
	8615A	Bipolar Forceps , jaws 45° upturned, for bipolar coagulation in larynx, for use with Bipolar High-Frequency Cord 847002V	23 cm	
	8615AS	Bipolar Forceps , jaws 45° upturned, with suction channel, for bipolar coagulation in larynx, for use with Bipolar High-Frequency Cord 847002V	23 cm	
	847002V	Bipolar High-Frequency Cord , for Valleylab coagulator, with two 2 mm cable sockets for Bipolar Suction Forceps 461010, 461015 and Bipolar Forceps 8615 A/AS	450 cm	
	8606E	Coagulation Suction Tube , with ergonomic handle, insulated, with connection for unipolar coagulation, with Cleaning Stylet 8606FM	23 cm	2.5 mm
	8606D	Coagulation Suction Tube , with ergonomic handle, insulated, with connection for unipolar coagulation, with Cleaning Stylet 8606FM	23 cm	2 mm
	8606F	Coagulation Suction Tube , with ergonomic handle, insulated, with connection for unipolar coagulation, with Cleaning Stylet 8606FM	23 cm	3.5 mm
	26006M-HN	Unipolar High-Frequency Cord , with 8 mm plug, for use with Valleylab	300 cm	
	26006M-D/10-HN	Unipolar High-Frequency Cord , with 8 mm plug, for use with Valleylab, sterile, for single use, package of 10	300 cm	

All product and company names are trademarks™ or registered® trademarks of their respective holders.


Clip Appliers (for use in Larynx)

Photo	Product Number	Description/Specifications	Length	Size
	8665L	LARYNGOFORCE® II Clip Forceps , jaws curved to left, with cleaning connector, for use with Ethicon LIGACLIP LT200-medium titanium clips only	22 cm	
	8665LL	LARYNGOFORCE® II Clip Forceps , curved 45° to left, for use with Ethicon LIGACLIP LT200-medium titanium clips only	22 cm	
	8665R	LARYNGOFORCE® II Clip Forceps , jaws curved to right, with cleaning connector, for use with Ethicon LIGACLIP LT200-medium titanium clips only	22 cm	
	8665RR	LARYNGOFORCE® II clip forceps , curved 45° to right, for use with Ethicon LIGACLIP LT200-medium titanium clips only	22 cm	

Clip Appliers & Hemostasis (for use in Diverticuloscopy & Esophagoscopy)

Photo	Product Number	Description/Specifications	Length	Size
	12067NL	Clip Forceps , jaws curved to left, with cleaning connector, for use with Ethicon Titanium Clips LT200	30 cm	
	12067NR	Clip Forceps , jaws curved to right, with cleaning connector, for use with Ethicon Titanium Clips LT200	30 cm	
	10390LN	Coagulation Suction Tube , insulated, with connector pin for unipolar coagulation, for use with Unipolar High-Frequency Cord 26002M	55 cm	4 mm
	10390BN	Coagulation Suction Tube , insulated, with connector pin for unipolar coagulation, for use with Unipolar High-Frequency Cord 26002M	35 cm	3 mm
	26006M	Unipolar High-Frequency Cord , with 8 mm plug, for use with Valleylab	300 cm	
	26006M-D/10-HN	Unipolar High-Frequency Cord , with 8 mm plug, for use with Valleylab, sterile, for single use, package of 10	300 cm	

Cleaning Brushes

Photo	Product Number	Description/Specifications
	27652-05-HN	Cleaning Brush , for cleaning jaws of CLICKLINE and TAKE-APART® instruments, for single use, package of 5

Notes

Notes

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.



KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34 | 78532 Tuttlingen/Germany
Postbox 230 | 78503 Tuttlingen/Germany
Phone: +49 7461 708-0 | Fax: +49 7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com

FOLLOW KARL STORZ

