


# Hysteroscopy – Enhancing Outcomes in Reproductive Medicine



Hysteroscopy may help improve pregnancy success rates through:

- Direct visualization of the inside of the uterine cavity
- Evaluation of the uterine lining
- Identification and specification of uterine abnormalities
- Treatment of uterine conditions that affect implantation and growth of the embryo

**STORZ**  
KARL STORZ – ENDOSKOPE

# Hysteroscopy – Enhancing Outcomes in Reproductive Medicine

Especially suitable for use in gynecological office settings and outpatient departments or IVF centers.

Hysteroscope	
CAMPO TROPHYscope®	
26008BAC	<b>CAMPO TROPHYscope®</b> , HOPKINS® Telescope 30°, diameter 2.9 mm
26152DA	<b>Continuous-Flow Examination Sheath</b> , size 3.7 mm
26152DB	<b>Continuous-Flow Operating Sheath</b> , size 4.4 mm

5 Fr. Operating Instruments	
26159UHW	<b>Biopsy and Grasping Forceps</b> , semirigid, double action jaws, 5 Fr., length 34 cm
26159SHW	<b>Scissors</b> , semirigid, pointed, single action jaws, 5 Fr., length 34 cm

Camera System	
TP101	<b>TELE PACK+</b>
TH110	<b>IMAGE1 S™ HX One-Chip FULL HD Camera Head</b>
495NTA	<b>Fiber Optic Light Cable</b> , diameter 2.5 mm, length 230 cm

Pump System	
UP220	<b>ENDOMAT® SELECT</b>
UP602	<b>Hysteroscopy Software</b>
031523-10	<b>Tubing Set</b> , irrigation

Mobile Stand	
UG803	<b>Subrack</b>
UG811	<b>Top Cover Monitor Holder</b>
UG854	<b>Stainless Steel Tube</b>
UG857	<b>Monitor Module</b>
UG860	<b>Drawer</b> , narrow



For more information scan the QR Code or email us at [gyn@karlstorz.com](mailto:gyn@karlstorz.com)

It is recommended to check the suitability of the product for the intended procedure prior to use.  
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.